## Drug Analysis Laboratory Quality Control – Balances 2012

Name/Model: Mettler MS201 S Manufacturer's Tolerance:

0 - 22g = +/- 0.0002g

Serial #: **B101100462** 

22-220g = +/-0.0003g

Location: Rm 362

	Jan	Feb	Mar	Apr	May	Jun	Jul	Aug	Sept	Oct	Nov	Dec
Day												
Chemist												
0.0020g*												
0.1000g												
1.0000g												
20.0000g												
50.0000g												
100.0000g												
200.0000g												
Sensitivity												
50.0000g and 10.0000g												
QC Review												
QA Review												

<sup>\*</sup> Balance used to make Standards.

Drug Analysis Laboratory Quality Control – Balances 2012

<sup>\*</sup> Balance used to make Standards.